



**NONESTERIFIED FATTY ACIDS (NEFA) REAGENTS
V514-0B**

CONTENTS

PRODUCT NO.

PACKAGING

NEFA REAGENTS	V514-0B	
ENZYME SAMPLE REAGENT (R1)	V514-01	3 x 8 mL
ENZYME COLOR REAGENT (R2)	V514-02	3 x 4 mL
REAGENT DILUENT	V514-03	2 x 25 mL

REAGENT PREPARATION

1. Reconstitute one vial of Enzyme Sample Reagent (R1) (Product No. V514-01) with 8 mL of Reagent Diluent (Product No. V514-03). Mix gently to prevent foaming until all material is in solution. Use within 24 hours.
2. Reconstitute one vial of Enzyme Color Reagent (R2) (Product No. V514-02) with 4 mL of Reagent Diluent (Product No. V514-03). Mix gently to prevent foaming until all material is in solution.
3. Store both reconstituted reagents at 2-8°C (refrigerated).

REAGENT STORAGE AND STABILITY

Store unopened and unreconstituted reagents at 2-8°C (refrigerated). When stored as directed, the reagents are stable until the expiration date stated on the label.

DIRECTIONS FOR USE

1. The Catachem NEFA test method requires two reagents (R1 and R2).
Use the reconstituted Enzyme Sample Reagent as R1 and the reconstituted Enzyme Color Reagent as R2.
2. After opening and reconstituting, use the R1 reagent within 24 hours.
Store refrigerated when not in use. R2 reagent after reconstitution are stable for 30 days when stored at 2-8°C (refrigerated). Cap tightly when not in use.
3. It is best to recalibrate the assay each day the assay is run.
4. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

**FOR TECHNICAL ASSISTANCE PLEASE CALL:
(203) 262-0330**

**OR EMAIL:
catachem@catacheminc.com
www.catacheminc.com**