



BILE ACIDS (CYCLING) REAGENT KIT
(with Calibrator and Bi-Level Controls) V562-0F

PRODUCT DESCRIPTION:

PRODUCT NO.

PACKAGING

BILE ACIDS (CYCLING) KIT

V562-0F

Sample Diluent (R1)

V562-10

1 x 6.5 mL

Activator Reagent (R2)

V562-11

1 x 1.5 mL

Bile Acids Calibrator

V404-23

1 x 1 mL

Bile Acids Control Level 1

V404-24

1 x 1 mL

Bile Acids Control Level 2

V404-25

1 x 1 mL

REAGENT PREPARATION:

The Bile Acids Sample Diluent Reagent (R1) and Activator Reagent (R2) are packaged ready for use. No preparation is required. The Bile Acids Calibrator, Control Level 1, and Control Level 2 are packaged ready for use. No preparation is required.

REAGENT STORAGE AND STABILITY:

Store the unopened Bile Acids Sample Diluent Reagent (R1), Activator Reagent (R2), Bile Acids Calibrator, Control Level 1 and Control Level 2 at 2-8°C. When stored as directed reagents, calibrator and controls are stable until the expiration date stated on the label. The reagents and controls have been tested to reflect shipping conditions and are stable for the lifespan of the product if frozen up to 5 times or reaching temperatures up to 40°C for up to one week.

DIRECTIONS FOR USE:

The Catachem, Inc. Bile Acids assay requires two reagents, R1 and R2. Use the Sample Diluent Reagent as R1 and the Activator Reagent R2. After opening, the Bile Acids Reagents are stable for at least 60 days at 2-8°C. Keep both reagents tightly capped while not in use. After opening, the Bile Acids Calibrator and Controls are stable for at least 30 days at 2-8°C. Keep tightly capped when not in use.
NOT FOR USE IN UNPROFESSIONAL SETTINGS.

FOR TECHNICAL ASSISTANCE PLEASE CONTACT:

Telephone: 203-262-0330

Or use the contact form at

www.catacheminc.com