



**IVDD, MHRA, CE**  
**ETHYLENE GLYCOL REAGENT KIT**  
**(including Calibrator and Bi-Level Controls)**  
**C504-0A**

[IN THE USA, FOR RESEARCH USE ONLY AND NOT FOR USE IN DIAGNOSTIC PROCEDURES]

**CONTENTS**

**PRODUCT NO.**

**PACKAGING**

**ETHYLENE GLYCOL KIT**

**C504-0A**

SAMPLE DILUENT (R1)

C504-01

2 x 10 mL

ACTIVATOR REAGENT (R2)

C504-02

1 x 8.5 mL

ETHYLENE GLYCOL CALIBRATOR

C504-10

1 x 3 mL

ETHYLENE GLYCOL CONTROL LEVEL I

C504-11

1 x 3 mL

ETHYLENE GLYCOL CONTROL LEVEL II

C504-12

1 x 3 mL

**REAGENT PREPARATION**

1. Add exactly 10 mL of distilled or deionized water to a bottle of Sample Diluent (R1), Product No. 504-01. Mix gently by inversion. The reconstituted Sample Diluent is now ready for use.
2. The Activator Reagent (R2), Product No. C504-02 and the Ethylene Glycol Calibrator, Control I and Control II are packaged ready for use. No preparation is required.

**REAGENT STORAGE AND STABILITY**

Store unopened Ethylene Glycol reagents, calibrator and controls at 2-8°C (refrigerated).  
When stored as directed, these materials are stable until the expiration date stated on the label.

**DIRECTIONS FOR USE**

1. The Catachem Ethylene Glycol Kit method test requires two reagents, R1 and R2.  
Use the reconstituted Ethylene Glycol Sample Diluent as R1 and the Ethylene Glycol Activator Reagent as R2. Once reconstituted, the Sample Diluent (R1) is stable on board the analyzer at 2-8°C (refrigerated) for 30 days if tightly capped when not in use.
2. After opening, the Ethylene Glycol Activator Reagent (R2), Calibrator, Control I and Control II are stable at 2-8°C (refrigerated) for 60 days if tightly capped when not in use.
3. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

**FOR TECHNICAL ASSISTANCE PLEASE CALL:**

**(203) 262-0330**

**OR EMAIL:**

**[catachem@catacheminc.com](mailto:catachem@catacheminc.com)**

**[www.catacheminc.com](http://www.catacheminc.com)**