

Intended Use


For **IN VITRO** quantitative determination of Lipase in serum using manual or automated applications.

Clinical Significance

The concentration of Lipase activity in blood is useful for the evaluation of pancreatic function. Lipase activity in serum is significantly elevated in cases of acute pancreatitis and obstruction of the pancreatic duct. Elevated values are also observed in intra-abdominal diseases, mumps and bacterial parotitis, but in these latter cases, values are generally lower than those seen in acute pancreatitis. Lower serum Lipase values have been observed in abscesses of the liver, acute hepatocellular damage, cirrhosis, cancer of the liver and bile duct. (1)

Method History

Human pancreatic Lipase catalyzes a two-step reaction: 1, 2 diglyceride is hydrolyzed to 2-monoglyceride and fatty acids. The 2-monoglyceride is further hydrolyzed by monoglyceride Lipase (MGLP) to glycerol and fatty acids. Glycerol is then quantitated using a first order rate GPO-Trinder reaction. Absorbance changes are monitored at 550 nm. The reaction scheme illustrates the reactions that take place in this method.



Reagent Content

When reconstituted according to the directions, the concentrations of the active ingredients in the reagent will be approximately as follows:

Enzyme Substrate

Each liter contains:	
1, 2 diglyceride	0.633g
Monoglyceride Lipase (microbial)	867 U
Glycerol kinase	1,333 U
Glycerol-3-phosphate oxidase (microbial)	40,000 U
TOOS	0.666g
ATP	0.400g
Peroxidase	1,333 U
Co-Lipase	40,000 U
Buffer	

Substrate Diluent

Each liter contains:	
Cholic acid	2.166g
Buffer	

Activator Reagent

Each liter contains:	
Deoxycholate	14.14 g
4-aminoantipyrine	1.20 g
Buffer	

Standard

Reconstitute each vial of standard/calibrator reagent with 3 ml of deionized water. Mix gently over a ten minute period. Use as directed.

Precautions

Handle the Catachem Lipase reagents following good clinical laboratory practice procedures. Avoid contact with skin and eyes. Should contact occur, wash affected area with plenty of cold water. **DO NOT PIPETTE REAGENTS BY MOUTH.**

Preparation of Working Reagents

Prepare the required number of vials of Catachem Lipase Working Enzyme Substrate Reagent by adding the contents of the 10 ml vial of Substrate Diluent to each Enzyme Substrate vial. Mix well and gently to prevent foaming, until completely in solution. Use the Activator Reagent as supplied. No preparation is required.

Reagent Storage and Stability

Store the Catachem Lipase reagents at 2-8°C. When stored as directed, these reagents are stable until the expiration date stated on the label. Store the Lipase Working Reagent (liquid) at 2-8°C. When prepared and stored as directed, the Working Reagent is stable for 21 days at 2°-8 °C.

Specimen Collection and Preparation

Test sera should be fresh, clear, and unhemolyzed. When blood is drawn, it should be processed as soon as possible and the serum should be isolated from the clot without delay. If the assays are not run immediately, the serum samples should be stored at 2-8°C or at -20°C.

Quality Control

To monitor the performance of the Working Reagents and the procedure used, we recommend the regular use of normal and abnormal control serum.

Interfering Substances

Several substances have been reported to alter the Lipase activity in serum. (5, 6) Glycerol concentrations greater than 100 mg/dL will produce erroneous results. Specimens with glycerol concentrations greater than 100 mg/dL should be diluted with physiological saline and reassayed.

Expected Values

The range of expected values in humans for this method is 7 - 59 u/L. These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the species under assay and for the area in which it is located.

Procedure

Important: Read entire procedure instructions before proceeding with assay.

Materials Required (Not Provided)

Spectrophotometer	
Cuvettes	1 cm light path
Timer	to time incubation time
Pipette	0.6 ml for reagent
Pipette	0.010 ml for sample

Materials Provided

Enzyme Substrate
Substrate Diluent
Activator Reagent
Standard:

Analytical Parameters

Wavelength	550 nm
Temperature	37°C
Pathlength	1 cm
Reaction Mode	Rate
Reaction Time	5 minutes
Reaction Volume (R1)	0.600 ml
Reaction Volume (R2)	0.200 ml
Sample Volume	0.010 ml
Total Volume	0.810 ml
Sample-to-reagent ratio	1:81

Assay Procedure

1. Prepare the required volume of Lipase Enzyme Substrate Working Reagent by following instructions for Working Reagent preparation.
2. Set spectrophotometer wavelength at 550 nm and zero the instrument with the cuvette containing water.
3. Pipette 0.6 ml of Working Reagent into each of two cuvettes marked: "Sample" and "Control".
4. Pipette 0.010 ml of control or sample into their respective cuvettes. Mix all cuvettes well.
5. Incubate cuvettes for 3.0 minutes at 37°C.
6. Pipette 0.200 ml of Activator Reagent into all cuvettes. Mix all cuvettes well and continuously monitor the change in absorbance for at least 3 minutes.
7. Read the "Control" and "Sample" absorbances.
8. Calculate the Lipase concentration (U/L) in the sample(s), as shown in calculations and results.

Calculations and Results

$$\text{Activity (U/L)} = \frac{\text{Sample (A2 - A1)}}{\text{Calibrator (A2 - A1)}} \times \text{Calibrator/standard (U/L)}$$

Example:	<u>A1</u>	<u>A2</u>
Sample	0.01	0.03
Calibrator	0.01	0.05

$$\begin{aligned} \text{Calibrator} &= 200 \text{ U/L} \\ &= \frac{0.03 - 0.01}{0.05 - 0.01} \times 200 \text{ U/L} \\ \text{Lipase (U/L)} &= 100 \text{ U/L} \end{aligned}$$

Units

One international unit (U/L) is defined as the amount of enzyme that catalyzes the conversion of one micromole of substrate per minute under the defined conditions

Method Performance Characteristics

Sensitivity: The sensitivity of this method is 0.008-0.0120 OD/U/L.

Linear Range: In this method there is no significant nonlinearity over the range of 0-600 U/L.

Precision Study

Lipase	SAMPLE 1	SAMPLE 2
N	15	15
U/L	18.80	72.60
SD	2.10	2.10
CV %	0.92	0.92

Correlation

A comparison of this method using an automated analyzer was carried out against a UV kinetic method. The performance of Catachem's Lipase method was found essentially equivalent to the reference method used. Linear regression analysis produced the following information:

$$\begin{aligned} N &= 58 \\ Y &= 1.844 x + 7.15 \\ r &= 0.982 \end{aligned}$$

References

1. Fundamentals of Clinical Chemistry. Edited by Norbert Tietz (1976) WB Saunders Co., Philadelphia.
2. Imamura S, Misaki H. A sensitivity method for assay of lipase activity by coupling with beta-oxidation enzymes of fatty acid. Selected Topics in Clinical Enzymology 2:73 (1984).
3. Imamura S, et al. A method using 1, 2-diglyceride as substrate for assay of lipase activity by coupling with beta-oxidation enzymes of fatty acid. Collection of summaries of lectures in the 126th general meeting of Kinki branch, analytical section, Japan Society of Clinical Chemistry pp 11-31 (1986).
4. Martin EW. Hazards of Medication. Alexander SF, Farage DJ, Hassan WE, Eds. Philadelphia PA and Toronto Canada. JB Lippincott Co (1971) pp 169-189.
5. Constantino NV, Kabat HF. Drug-induced values-revised 1973. Am J Hosp Pharm 30:24-71 (1973).
6. Young DS, Pestaner LC, Gibberman V. Effects of drugs on clinical laboratory tests. Clin Chem 21 No 5 (1975).

REV: BK02200211dt0513