

Assay Procedure

1. Label cuvettes or appropriate test tubes as: a) Calibrator Blank (CAL-BLK), b) Calibrator (CAL), c) Control 1 Blank (C-1BLK), d) Control 1 (C1), e) Control 2 Blank (C-2BLK), f) Control 2 (C2), g) Sample Blank (SAMP BLK), h) Sample (SAMP).
2. Pipette the reagent and sample volumes into the cuvettes or test tubes as shown in table below. Pipette Catachem NEFA Enzyme Reagent I (R-1) first, followed by the sample.
3. Incubate for 5 minutes.
4. Immediately, after the incubation period is over, add Catachem NEFA Enzyme Reagent II (R-2) to the assay samples and water to the blank samples. Mix all cuvettes.
5. Set a timer for exactly 4 minutes.
6. At the end of the 4 minutes read all cuvettes at 550nm. Record all absorbencies.

	CAL BLK	CAL	C-1 BLK	C-1	C-2 BLK	C-2	SAMP BLK	SAMP
	ml	ml	ml	ml	ml	ml	ml	ml
RGT 1	0.5	0.5	0.5	0.5	0.5	0.5	0.5	0.5
SAMP	0.050	0.050	0.050	0.05	0.050	0.05	0.050	0.050
INCUBATE FOR 5 MINUTES								
RG2 2	0.00	0.10	0.00	0.10	0.00	0.10	0.00	0.10
H2O	0.10	0.0	0.10	0.00	0.10	0.00	0.10	0.00
MIX AND INCUBATE FOR 4 MINUTES								
MIX AND READ ALL CUVETTES								

Calculations and Results

$$\text{NEFA mmol/L} = \frac{\Delta\text{- Abs. Samp.}}{\Delta\text{- Abs. Cal.}} \times \text{Cal. Value (mmol/L)}$$

Example:

	<u>Samp. Abs.</u>	<u>Blk. Abs.</u>	<u>Δ-Abs.</u>
Sample	0.400	0.300	0.100
Calibrator	0.150	0.030	0.120

Calibrator assigned value = 0.5 mmol/L

$$\begin{aligned} \text{Sample NEFA mmol/L} &= \frac{0.100}{0.120} \times 0.5 \text{ mmol/L} \\ &= 0.42 \text{ mmol/L} \end{aligned}$$

Quality Control

To monitor the quality performance of the procedure used, Catachem NEFA Control Level I and Control Level II should be included in the assay procedure each time the assay is run.

Procedure Limitations

Samples with NEFA values greater than 2.5mmol/L should be diluted 1:2 with physiological saline and reassayed. Multiply results obtained by 2 to adjust for the sample dilution.

Interfering Substances

The following substances have no significant effect on the accuracy of this NEFA procedure at the concentrations stated.

- Hemoglobin ≤ 200 mg/dL
- Bilirubin 10.0 mg/dL
- Ascorbic Acid ≤ 20 mg/dL

Other substances and certain drugs are also known to influence the NEFA values (6).

Method Performance Characteristics

Sensitivity: Using a path length of 1 cm, a Δ-absorbance of 0.1-0.20 per mmol/L should be obtained.

Linearity: This procedure is linear over the range of 0-2.5 mmol/L.

Precision: Precision data was obtained using five levels of protein based controls and following the NCCLS EP5-T2 procedure (5). The following results were observed:

Precision

NEFA Mean mmol/L	Within-Run Precision		Total Precision	
	SD mmol/L	CV %	SD mmol/L	CV %
0.169	0.0060	3.550	0.0074	4.379
0.443	0.0087	1.953	0.0100	2.257
1.02	0.0149	1.461	0.0162	1.588
1.516	0.0163	1.078	0.0172	1.134

ACCURACY

Using an automated analyzer, correlation studies were carried out between this Catachem NEFA procedure (Y) and a commercially available NEFA test kit as reference (X). Serum samples were assayed and the results compared by the least squares regression. The following statistics were observed:

N	=	44
Range	=	0.1-2.3 mmol/L
Mean Y	=	0.5295 mmol/L
Mean X	=	0.5841 mmol/L
Y	=	0.9137x - 0.0042
r	=	0.9949
Sy.x	=	0.0541

Bibliography

1. Fundamentals of Clinical Chemistry. Edited by Norbert Tietz 2nd. Ed. Philadelphia: W.B Saunders; 1976.
2. Itaya, K and Ui, M. Colorimetric determination of free fatty acids in biological fluids. J. Lipids Research; 6, 16 (1965).
3. Novak, M. Colorimetric ultramicro method for the determination of fatty acids. J. Lipid Research; 6, 431 (1965).
4. Trout, D.L., Estes, E.H. and Friedberg, S.J. J. Lipid Research; 1,199 (1960).
5. Evaluation of Precision Performance of Clinical Chemistry Devices. Second Edition. NCCLS Document EP5-T2. Vol. 12, No. 4.
6. Young D.S., Pestamer L.C., Gibberman V. Effects of drugs on clinical laboratory tests. Clin. Chem. 21, No. 5(1975).

REV: 051215DT