



### Intended Use


For **In Vitro Diagnostic** use in the quantitative determination of Pyruvate in serum or plasma.

### Clinical Significance (1)

Pyruvate is one of main intermediates in the metabolic pathway. It can be converted into energy through acetyl-CoA, to a carbohydrate through gluconeogenesis, to the amino acid alanine and to ethanol. Unusual pyruvate levels can be found in various metabolic disorders and also liver diseases.

### Method Principle

Pyruvate, by the action of Pyruvate Oxidase converts pyruvate in the presence of other components to acetyl phosphate with the release of hydrogen peroxide. The hydrogen peroxide thus produced is quantitatively determined by coupling 4-aminoantipyrine with N-ethyl-N-(2-hydroxy-3-sulphopropyl)-m-toluidine (TOOS) (3-4), where a quinonemine dye with maximum absorption at between 540-550 nm is produced (2). The following reaction scheme illustrates the reactions that occur in this method:



### PYRUVATE REAGENT

Each liter contains:

Buffer	
4-Aminoantipyrine (4AA)	0.15 mol
TOOS	1.2 mmol
Pyruvate Oxidase	≥ 200 Units
Peroxidase	≥ 5000 Units
Stabilizers, co-factors and nonreactive ingredients	

### Precautions

Handle this reagent using good laboratory practice. **DO NOT PIPETTE REAGENT BY MOUTH.** Avoid contact with skin and eyes. If contact occurs, wash affected area with plenty of cold water. Clean spills immediately.

### Reagent Storage and Stability

Store the Catachem Pyruvate reagent at 2-8°C. When stored as directed the reagent is stable until expiration date stated on the label.

### Working Reagent Preparation

The Catachem Pyruvate reagent is packaged in a powder format which is reconstituted with deionized/distilled water to form a single working reagent. Once opened and reconstituted, the Catachem Pyruvate Reagent is stable for at least 14 days at 2-8°C if stored capped.

### Reagent Indications of Deterioration

- Turbidity
- Absorbance > 0.5 OD, 1 cm light path, 540/550nm
- Quality control values out of assigned ranges

If these reagent characteristics are observed contact Catachem technical service.

### Specimen Collection and Stability (2, 3)

To maintain sample integrity and avoid changes in Pyruvate concentrations, care should be taken when collecting sample specimens.

Venous specimens should be collected without the use of a tourniquet or immediately after a tourniquet has been applied. Literature recommends sample collection in 6% perchloric acid in a chilled collection tube. (6)

Specimens should be collected, processed and analyzed promptly.

### Procedure

These instructions are outlined for performing the Pyruvate assay by manual or automated procedures. Read the original instrument manufacturer's instructions and procedures before performing this as an automated pyruvate procedure.

### Materials Provided

Catachem Pyruvate Reagent

### Materials required but not Provided

- Analyzer (spectrophotometer) with 540 - 550nm wavelength
- Calibrator material with assigned pyruvate value
- Quality control material with assigned pyruvate values
- Deionized/distilled water

### Calibration

Catachem Pyruvate Calibrator with assigned pyruvate value is recommended.

### Calibration Schedule

Calibration should be performed when this method is run manually or implemented on an automated analyzer for the first time. Recalibration is required after changes of reagent lot number, major instrument service, and when quality control values are out of the indicated range.

### Calibration Procedure

Instructions for calibrating any automated analyzer are provided by the specific instrument manufacturer. Read the entire recommended calibration procedure before proceeding with the instrument calibration.

### Quality Control

To monitor the quality performance of the procedure, Catachem Pyruvate Control Level I and Level II with assigned pyruvate ranges should be included in the assay procedure each time the assay is performed.

### Directions for Use

The Catachem Pyruvate method requires one reagent.

### Procedure Limitations

Samples with Pyruvate values greater than 700 µmol/L should be diluted 1:2 with physiological saline and reassayed. Multiply results obtained by 2 to adjust for the sample dilution.

### Interfering Substances

The following substances have no significant effect on the accuracy of this Pyruvate procedure at the concentrations stated.

- Hemoglobin ≤ 200 mg/dL
- Triglycerides ≤ 1000 mg/dL
- Bilirubin ≤ 2.2 mg/dL
- Lactate ≤ 200 µmol/L

Other substances and certain drugs may also influence the pyruvate values. (4)

**Normal values (5)**

In humans – 34-102 µmol/L (0.3-0.9 mg/dL)  
 In other animals, values vary depending on species

**Analytical Parameters (Manual method)**

Wavelength	540nm
Temperature	37° C
Path length	1 cm
Reaction Mode	End Point
Reaction Time	8 minutes
Reaction Volume	0.5 mL
Sample Volume	0.025 mL
Total Volume	0.525 mL
Sample-to-reagent ratio	1:20

**Assay Procedure**

1. With deionized/distilled water reconstitute the required volume of Pyruvate Working Reagent bottles needed dependent on the number of samples to be assayed.
2. Set spectrophotometer temperature at 37°C, wavelength at 540nm.
3. Pipette 0.5 mL of Pyruvate reagent into each cuvette marked "Blank", "Test", "Calibrator" and "Control".
4. Pipette 0.025 mL of Blank (water), Test sample, Calibrator or Control into their respective cuvettes.
5. Mix all cuvettes and incubate all cuvettes for 5 minutes.
6. Read absorbance (A) of each cuvette marked "Blank", "Test", "Calibrator" and "Control" at the set wavelength of 540nm.

Pyruvate Procedure Scheme				
	RGT BLANK	CAL	CONTROL	TEST
Pyruvate Rgt.	0.5 mL	0.5 mL	0.5 mL	0.5 mL
Water (Blank)	0.025 mL	-	-	-
Sample	-	0.025 mL	0.025 mL	0.025 mL
Mix	All cuvettes			
Incubate	All cuvettes at 37°C for 5 minutes			
Calculate	Pyruvate µmol/L			

**Calculations and Results:**

$$\text{Pyruvate } (\mu\text{mol/L}) = \frac{(\text{A}) \text{ Sample} - (\text{A}) \text{ Blank}}{(\text{A}) \text{ Calibrator} - (\text{A}) \text{ Blank}} \times \text{Calibrator } (\mu\text{mol/L})$$

CATACHEM PYRUVATE REAGENT	
<b>General BECKMAN AU TEST PARAMETERS</b>	
Predilution Rate	1
Sample	Serum
S Volume	7.5
S Diluent	0
Repeat S Vol	*
R1 Vol	150
R1 Diluent	0
R2 Vol	-
R2 Diluent	0
Wave Length Primary	540 nm
Wave Length Secondary	700 nm
Method Type/Reaction Slope	END +
Read Point Start 1	0
Read Point End 1	27
Read Point Start 2	-

Read Point End 2	-
Linearity Limit	*
Reaction Min OD	-2.0
Reaction Max OD	3.0
Reagent OD Range 1-L	-2
Reagent OD Range 1-H	2.5
Reagent OD Range 2-L	-2
Reagent OD Range 2-H	2.5
Dynamic Range Low	0
Dynamic Range High	700
Normal L	*
Normal H	*
On Board Stability (Days #)	14**
Cal Factor Low	0
Cal Factor High	99999
Calibration Type	AB
Formula	Y=AX+B
Calibration Replicates	
Cal 1	*
Cal 2	
Cal Frequency Days	*
Units	µmol/L

\*User defined \*\* If Calibrated before each patient run

**Method Performance Characteristics**

**Linearity:** This procedure is linear over the range of 0-700 µmol/L

**Bibliography**

1. Tietz NW. ed. Fundamentals of Clinical Chemistry. 2<sup>ND</sup> ed. Philadelphia. WB Saunders, 1976: 936-940.
2. Trinder P, Ann Clin Biochem 6, 24 (1969).
3. Kleine TO. Nervensysteme. In: Greiling H, Greesner AM, eds. Lehrbuch der Klinischem Chemie und Pathobiochemie. Stuttgart: Schattauer: 1987: 839-893.
4. Young D.S., Pestamer L.C., Gibberman V. Effects of drugs on clinical laboratory tests. Clin. Chem. 21, No. 5(1975).
5. Hansen J.L. & Freier E.F. Clin Chem. 24/3, 4755-479 (1978)
6. Mayo Medical Laboratories -CPT code 84210 - Pyruvate

Rev -20220728sddmt