



<u>PRODUCT DESCRIPTION</u>	<u>PRODUCT NO.</u>	<u>PACKAGING</u>
METHANOL KIT	C604-0A	
METHANOL REAGENT	C604-01	2 x 10 mL [Lyo]
METHANOL REAGENT DILUENT	C604-02	2 x 10 mL
METHANOL CALIBRATOR	C604-10	1 x 3 mL
METHANOL CONTROL LEVEL I	C604-11	1 x 3 mL
METHANOL CONTROL LEVEL II	C604-12	1 x 3 mL

REAGENT PREPARATION

1. Prepare Catachem Methanol Working Reagent by pouring the entire contents (10 mL) of Reagent Diluent (Product No. C604-02) into the Methanol Reagent bottle (Product No. C604-01). Mix gently by inversion for 3-5 minutes. The Methanol Calibrator, Control Level 1 and Control Level 2 are ready to use and require no preparation.
2. After preparation, the Methanol Working Reagent is stable at 2-8°C for 30 days. After opening, the Calibrator, Control I and Control II are stable at 2-8°C for 60 days.

REAGENT STORAGE AND STABILITY

Store the unopened Methanol Reagents at 2-8°C. When stored as directed, these reagents are stable until the expiration date stated on the label. Store the unopened Methanol Calibrator and Controls at 2-8°C. When stored as directed, the Methanol Calibrator and Controls are stable until the expiration date stated on the label.

DIRECTIONS FOR USE

The Catachem Methanol method requires one working reagent. The Working Reagent is stable aboard the analyzer at 2-8°C for 30 days. Cap working reagent bottle when not in use.

FOR TECHNICAL ASSISTANCE PLEASE CALL:

**203-262-0330
OR EMAIL**

**Catachem@catacheminc.com
www.catacheminc.com**