



CARBON DIOXIDE (CO₂) REAGENT
C206-02

Contents	Product No.	Package
CARBON DIOXIDE (CO₂) REAGENT (Cycling)	C206-02	6 x 50 mL

REAGENT PREPARATION

Catachem Carbon Dioxide (cycling) is a liquid, ready to use reagent.
No preparation is required.

REAGENT STORAGE AND STABILITY

Store unopened reagents at 2-8°C.
When stored as directed, the reagent is stable until the expiration date stated on the label.

NOT FOR USE IN UNPROFESSIONAL SETTINGS

FOR TECHNICAL ASSISTANCE:
Email: catachem@catacheminc.com
Contact Form: www.catacheminc.com
Call: 203-262-0330



CARBON DIOXIDE (CO₂) REAGENT C206-02 MANUAL/AUTOMATED APPLICATION

Intended Use

For the quantitative determination of Carbon Dioxide (Bicarbonate) in serum or plasma using a manual application.

Clinical Significance



Determinations of Bicarbonate are primarily used for diagnosing respiratory acidosis and metabolic alkalosis, as well as for monitoring the causes and treatments. ⁽¹⁾

Method History

In 1921 Van Slyke and Stadie ⁽²⁾ published a method for the determination of gases in blood. With availability of phosphoenolpyruvate carboxylase (PEPC), an enzyme which utilizes Bicarbonate as a substrate, simpler and specific methods have been reported. ^(3,4) The Catachem Carbon Dioxide (Bicarbonate) procedure for manual or automated applications utilizes a PEPC extracted from a specific microorganism and an enzyme cycling system to prevent contamination by CO₂ from the air.

Method Principle

In the presence of bicarbonate, PEPC catalyzes the carboxylation of PEP to form oxaloacetate. Oxaloacetate is then quantified via the MDH-Acetyl NADH reaction with concomitant oxidation of Acetyl NADH to Acetyl NAD. The decrease of the Acetyl NADH absorbance is monitored at 380nm and it is directly proportional to the concentration of bicarbonate in the sample. The reaction scheme illustrates the reactions that occur in this method.



Reagent

The reagent comes as a stable liquid. The concentrations of the active ingredients in the reagents are approximately as follows:

Carbon Dioxide Reagent

Each liter contains:

Buffer	
Phosphoenolpyruvate	3.5 mmol
Acetyl NADH	0.6 mmol
PEPC (Microbial)	≥ 600 units
MDH (Microbial)	≥ 750 units
Non-reactive ingredients	

Precautions

Avoid contact of reagents with skin and eyes. Should contact occur, wash affected area with plenty of water. **DO NOT PIPETTE REAGENTS BY MOUTH.**

Reagent is provided ready to use.

Reagent Storage and Stability

Store Catachem CO₂ Reagent at 2-8°C. When stored unopened as directed, the reagent is stable until the expiration date stated on the

label. Once opened the reagent is stable for at least **60 days** or until the baseline OD at 380 nm is less than < 0.5 OD units. After opening when not in use, it is recommended to store the reagent in a **tightly capped** container. The Catachem CO₂ Reagent has been tested to reflect shipping conditions and is stable for the lifespan of the product if frozen up to 5 times or upon reaching temperatures up to 40°C for up to one week.

Specimen Collection and Preparation

Serum samples that have been stored at room temperature and exposed to ambient air are not recommended since with time, a significant amount of CO₂ will be lost. Samples stored at 2-8°C and tightly sealed are stable for approximately one hour. ⁽⁴⁾

Quality Control

To monitor the performance of the Working Reagents and the procedure used, we recommend the regular use of a normal and abnormal control serum.

Interfering Substances

Several substances have been reported to interfere with the CO₂ method. ⁽³⁾ A summary of interference from drugs on clinical laboratory procedures may be found by consulting D.S. Young, et al. ⁽⁵⁾

Expected Values

The range of expected values for humans determined for this method is 22 mEq/L to 28mEq/L. Ranges for animals vary by species. All values are considered suggested guidelines. It is recommended that each laboratory establish the normal range for the species under test and for the area in which it is located. ⁽⁴⁾

Procedure

Important: Read entire procedural instructions before proceeding with assay.

Materials Required (Not Provided)

Spectrophotometer	
Match cuvettes	1 cm light path
Timer	to time incubation time
Pipette	1.0 ml for reagent
Pipette	0.01 ml for sample
Calibrator	Catachem Catacal or equivalent calibrator

Material Provided

Catachem CO₂ Reagent

Analytical Parameters

Wavelength	380 nm
Pathlength	1 cm
Reaction mode	endpoint
Reaction time	5.0 minute
Reagent volume	1.0 ml
Sample volume	0.010 ml
Total volume	1.01 ml
Sample to reagent ratio	1:101



CARBON DIOXIDE (CO₂) REAGENT
C206-02
MANUAL/AUTOMATED APPLICATION

Assay Procedures

1. Pipette 1.0 ml of CO₂ reagent into each of three cuvettes marked "Calibrator", "Sample", and "Blank".
2. Pipette 0.01 ml of Calibrator or Sample into their respective cuvettes. Use 0.01 ml of distilled water for the Blank. Mix all cuvettes well.
3. Incubate all cuvettes for 5.0 minutes at 37°C.
4. Set spectrophotometer wavelength at 380nm and zero the instrument with the Blank.
5. Read the "Calibrator" and "Sample" absorbencies.
6. Calculate the CO₂ concentration (meq/L) in the Sample(s) as shown in calculations and results.

Calculations and Results (note Δ means "change in".)

$$\text{Bicarbonate (meq/L)} = \frac{\text{Sample } \Delta \text{ OD}}{\text{Calibrator } \Delta \text{ OD}} \times \text{Calibrator (meq/L)}$$

Example:

Sample Δ OD	0.300
Calibrator Δ OD	0.200
Calibrator	= 20 meq/L

$$\text{Bicarbonate (meq/L)} = \frac{0.30}{0.20} \times 20 \text{ meq/L}$$

$$= 30 \text{ meq/L}$$

Method Performance Characteristics

Sensitivity: The sensitivity of this method is 0.023-0.027 absorbance units/meq/L.

Linear Range: In this method there is no significant nonlinearity over the range of 0-45 meq/L.

Precision: Within-run and day-to-day precision is summarized below:

Precision Study

CO ₂	Within-Run		Total Precision	
	SD	CV	SD	CV
Meq/L	Meq/L	%	Meq/L	%
16.0	0.11	0.71	0.15	0.96
23.5	0.27	1.15	0.50	2.14
29.8	0.21	0.70	0.49	1.67

Correlation

A comparison of this method using an automated analyzer and a reference method based upon the PEPC reaction resulted in the following regression statistics.

Range	=	9-36 meq/L
N	=	122
Y	=	1.05x + 0.38
r	=	0.977
Sy.x	=	1.23

References

1. Fundamentals of Clinical Chemistry. Edited by Norbert Teitz. WB Saunders, Philadelphia (1976).
2. Van Slyke DD and Stadie WC. J. Biol. Chem. 49.1 (1921).
3. Menson RC, Narayanswamy V, Bussian RW and Adams TH. Clin Chem 20, 872 (1972).
4. Tietz NW. Textbook of Clin Chem. WB Saunders Co. pp1188, 1209-1210 (1986).
5. Young GS, Pestaner LC and Gibberman V. Effect of drugs on clinical laboratory tests. Clin Chem 21 1D (1975).