



*Vet Spec*  
**ALBUMIN REAGENT**  
**V244-12; C244-50, C244-01**

<u>CONTENTS</u>	<u>PRODUCT NO.</u>	<u>PACKAGING</u>
ALBUMIN REAGENT	V244-12	6 x 13 mL
ALBUMIN REAGENT	C244-50	1 x 50 mL
ALBUMIN REAGENT	C244-01	1 x 500 mL

### REAGENT PREPARATION

This reagent is packaged ready for use. No preparation is required.

### REAGENT STORAGE AND STABILITY

Store unopened reagent at 15 - 30°C (room temperature) or at 2 - 8 °C (refrigerated). When stored as directed, the reagent is stable until the expiration date stated on the label.

### DIRECTIONS FOR USE

The Catachem Albumin Reagent, after opening, is stable until the expiration date stated on the label. Opened reagent may also be stored on board an automated analyzer at 2 - 8 °C (refrigerated) or at 8 - 15 °C. Cap reagent when not in use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

**FOR TECHNICAL ASSISTANCE PLEASE CALL:**

**(203) 262-0330**

**OR EMAIL:**

**[catachem@catacheminc.com](mailto:catachem@catacheminc.com)**

**[www.catacheminc.com](http://www.catacheminc.com)**