



Vet Spec
GLUCOSE REAGENT
C124-12, C124-06, C124-07

CONTENTS

PRODUCT NO.

PACKAGING

GLUCOSE REAGENT

C124-12

6 x 13 mL

GLUCOSE REAGENT

C124-06

6 x 25 mL

GLUCOSE REAGENT

C124-07

6 x 50 mL

REAGENT PREPARATION

This reagent is packaged ready for use. No preparation is required.

REAGENT STORAGE AND STABILITY

Store unopened reagent at 2-8°C (refrigerated). When stored as directed, the reagent is stable until the expiration date stated on the label.

DIRECTIONS FOR USE

The Catachem Glucose Reagent, after opening, is stable on board most analyzers for at least 60 days at 2-8°C (refrigerated). Cap tightly when not in use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

FOR TECHNICAL ASSISTANCE PLEASE CALL:

(203) 262-0330

OR EMAIL:

catachem@catacheminc.com

www.catacheminc.com