



**CONTENTS**

**PRODUCT NO.**

**PACKAGE**

<b>AST KIT – A</b>	<b>C154-0A</b>	
<b>AST SAMPLE DILUENT</b>	<b>C154-03</b>	<b>6 x 22.5ml</b>
<b>AST ACTIVATOR</b>	<b>C154-04</b>	<b>6 x 2.5ml</b>
<b>AST KIT – B</b>	<b>C154-0B</b>	
<b>AST SAMPLE DILUENT</b>	<b>C154-05</b>	<b>6 x 45.0ml</b>
<b>AST ACTIVATOR</b>	<b>C154-06</b>	<b>6 x 5.0ml</b>

**REAGENT PREPARATION**

The Catachem AST reagents are packaged ready for use. No preparation of either reagent is required.

**REAGENT STORAGE AND STABILITY**

Store the unopened reagent kit at 2 - 8°C (refrigerated). When stored as directed, both unopened reagents are stable until the expiration date stated on the label.

**DIRECTIONS FOR USE**

The Catachem AST test method requires two reagents: the Diluent as R1 and the Activator as R2. After opening, these reagents are stable at 2 - 8°C for at least 30 days. Cap when not in use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

**FOR TECHNICAL ASSISTANCE PLEASE CALL:**

**203-262-0330**

**OR EMAIL:**

**[catachem@catacheminc.com](mailto:catachem@catacheminc.com)**

**[www.catacheminc.com](http://www.catacheminc.com)**