



**ALANINE AMINOTRANSFERASE (ALT) REAGENT
C164-0A, C164-0B**

<u>CONTENTS</u>	<u>PRODUCT NO.</u>	<u>PACKAGE</u>
ALT KIT – A	C164-0A	
ALT SAMPLE DILUENT	C164-03	6 x 22.5ml
ALT ACTIVATOR	C164-04	6 x 2.5ml
ALT KIT – B	C164-0B	
ALT SAMPLE DILUENT	C164-05	6 x 45.0ml
ALT ACTIVATOR	C164-06	6 x 5.0ml

REAGENT PREPARATION

The Catachem ALT reagents are packaged ready for use. No preparation of either reagent is required.

REAGENT STORAGE AND STABILITY

Store the unopened reagent kit at 2 - 8°C (refridgerated). When stored as directed, both unopened reagents are stable until the expiration date stated on the label.

DIRECTIONS FOR USE

The Catachem ALT test method requires two reagents: the Diluent as R1 and the Activator as R2. After opening, these reagents are stable at 2 - 8°C for at least 30 days. Cap when not in use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

FOR TECHNICAL ASSISTANCE PLEASE CALL:

203-262-0330

Or use the contact form at

www.catacheminc.com