



**ALANINE AMINOTRANSFERASE (ALT) REAGENT
C164-0A, C164-0B**

| <u>CONTENTS</u> | <u>PRODUCT NO.</u> | <u>PACKAGE</u> |
|--------------------|--------------------|----------------|
| ALT KIT – A | C164-0A | |
| ALT SAMPLE DILUENT | C164-03 | 6 x 22.5ml |
| ALT ACTIVATOR | C164-04 | 6 x 2.5ml |
| ALT KIT – B | C164-0B | |
| ALT SAMPLE DILUENT | C164-05 | 6 x 45.0ml |
| ALT ACTIVATOR | C164-06 | 6 x 5.0ml |

REAGENT PREPARATION

The Catachem ALT reagents are packaged ready for use. No preparation of either reagent is required.

REAGENT STORAGE AND STABILITY

Store the unopened reagent kit at 2 - 8°C (refridgerated). When stored as directed, both unopened reagents are stable until the expiration date stated on the label.

DIRECTIONS FOR USE

The Catachem ALT test method requires two reagents: the Diluent as R1 and the Activator as R2. After opening, these reagents are stable at 2 - 8°C for at least 30 days. Cap when not in use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

FOR TECHNICAL ASSISTANCE PLEASE CALL:

203-262-0330

Or use the contact form at

www.catacheminc.com