



CONTENTS

PRODUCT

PACKAGING

CREATININE REAGENT KIT	C324-0A	
DILUENT REAGENT (R1)	C324-01	1 x 250 mL
COLOR REAGENT(R2)	C324-02	1 x 250 mL

REAGENT PREPARATION

These reagents are packaged ready for use. No preparation is required.

REAGENT STORAGE AND STABILITY

Store the unopened Diluent Reagent (R1) and unopened Color Reagent (R2) at 15-30°C (room temperature) or at 2-8°C (refrigerated). When stored as directed, both reagents are stable until the expiration date stated on the label.

DIRECTIONS FOR USE

The Catachem Vet Spec Creatinine test method requires two reagents: Diluent Reagent as R1 and Color Reagent as R2. After opening, these reagents are stable until the expiration date stated on the label. Opened reagents may also be stored on board an automated analyzer at 2-8°C (refrigerated) or at 8-15°C. Cap reagents when not is use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

FOR TECHNICAL ASSISTANCE PLEASE CALL:

(203) 262-0330

OR EMAIL:

catachem@catacheminc.com

www.catacheminc.com