



*Vet Spec*  
**MAGNESIUM REAGENT**  
**V355-12, C355-01**

| <u>CONTENTS</u> | <u>PRODUCT NO.</u> | <u>PACKAGING</u> |
|-----------------|--------------------|------------------|
| MAGNESIUM       | V355-12            | 6 x 13 mL        |
| MAGNESIUM       | C355-01            | 1 x 125 mL       |

### REAGENT PREPARATION

This reagent is packaged ready for use. No preparation is required.

### REAGENT STORAGE AND STABILITY

Store unopened reagent at 15 - 30°C (room temperature) or at 2 - 8°C (refrigerated). When stored as directed, the reagent is stable until the expiration date stated on the label.

NOTE. Magnesium reagent deteriorates quickly if exposed to air.

### DIRECTIONS FOR USE

The Catachem Vet Spec Magnesium Reagent, after opening, is stable for up to 30 days if capped when not in active use. Opened reagent may also be stored on board an automated analyzer at 2 - 8°C (refrigerated) or at 8 - 15°C for up to 30 days if capped when not in use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

**FOR TECHNICAL ASSISTANCE PLEASE CALL:**

**(203) 262-0330**

**OR EMAIL:**

**[catachem@catacheminc.com](mailto:catachem@catacheminc.com)**

**[www.catacheminc.com](http://www.catacheminc.com)**