



| <u>PRODUCT DESCRIPTION:</u> | <u>PRODUCT NO.</u> | <u>PACKAGING</u> |
|-----------------------------|--------------------|------------------|
| Fructosamine Kit | C412-0A | 3 x 50 mL |
| Fructosamine Reagent | C414-02 | 1 X 50 mL |
| Fructosamine Reagent | V412-12 | 3 x 13 mL |

REAGENT PREPARATION:

This reagent is packaged in ready-to-use form. Mix gently before use.

REAGENT STORAGE AND STABILITY:

Store unopened reagent at 2-8°C (refrigerated). When stored as directed, the reagent is stable until the expiration date stated on the label. The reagent has been tested to reflect shipping conditions and is stable for the lifespan of the product if frozen up to 5 times or reaching temperatures up to 40°C for up to one week.

DIRECTIONS FOR USE:

The Catachem, Inc. Fructosamine Reagent, after opening, is stable on board the analyzer for a minimum of 6 months at 2-8°C (refrigerated) and capped when not in use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

FOR TECHNICAL ASSISTANCE PLEASE CONTACT:

Telephone: 203-262-0330

Or use the contact form at

www.catacheminc.com