



**PRODUCT DESCRIPTION:**

**PRODUCT NO.**

**PACKAGING**

Bromide Sample Diluent	C424-01	1 x 30 mL
Bromide Color Reagent	C424-02	1 x 15 mL

**BROMIDE KIT**

**C424-0B**

Bromide Sample Diluent	C424-01	1 X 30 mL
Bromide Color Reagent	C424-02	1 X 15 mL

**BROMIDE KIT**

**C422-0A**

Bromide Sample Diluent	C422-01	3 x 60 mL
Bromide Color Reagent	C422-02	3 x 30 mL

**REAGENT PREPARATION:**

Both reagents are packaged in ready-to-use form. No preparation is required.

**REAGENT STORAGE AND STABILITY:**

Store both Reagents at 2-8°C. When stored as directed, both reagents are stable until the expiration date stated on the label. The reagents have been tested to reflect shipping conditions and are stable for the lifespan of the product if frozen up to 5 times or reaching temperatures up to 40°C for up to one week.

**DIRECTIONS FOR USE:**

The Bromide method requires two reagents, R1 and R2. Use the Bromide Sample Diluent as R1 and the Bromide Color Reagent as R2. Once opened, the Sample Diluent (R1) is stable at 2-8°C until the expiration date stated on the label. Once opened, the Color Reagent (R2) is stable for 60 days at 2-8°C. Both reagents should be tightly capped when not in use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

**FOR TECHNICAL ASSISTANCE PLEASE CONTACT:**

**Telephone: 203-262-0330**

**Or use the contact form at**

**[www.catacheminc.com](http://www.catacheminc.com)**