



**GLUTAMATE DEHYDROGENASE REAGENT
POWDER SAMPLE
C552-0A**

CONTENTS

PRODUCT NO.

PACKAGE

| | | |
|-------------------------|----------------|-------------|
| GLDH REAGENT KIT | C552-0A | |
| GLDH REAGENT DILUENT | C552-01 | 3 x 25 mL |
| GLDH REAGENT POWDER | C552-02 | 3 x 67.5 mg |

REAGENT PREPARATION

To prepare the Catachem GLDH powder reagent empty the entire contents of one bottle of powder into one bottle of diluent. Swirl to combine. Once activated the GLDH powder reagent is stable for at least 5 days when stored at 2-8°C and capped tightly when not in use.

DIRECTIONS FOR USE

Unopened and stored at 2-8°C, the GLDH reagent diluent and powder are stable until the expiry date stated on the bottle labels. Cap tightly when not in use.

NOT FOR USE IN UNPROFESSIONAL SETTINGS.

FOR TECHNICAL ASSISTANCE PLEASE CALL:

(203) 262-0330

catachem@catacheminc.com

www.catacheminc.com