



CHOLESTEROL REAGENT V104-12, C104-01,02, C103-01

CONTENTS

PRODUCT NO.

PACKAGING

CHOLESTEROL REAGENT	V104-12	6 x 13 mL
CHOLESTEROL REAGENT	C104-01	6 x 25 mL
CHOLESTEROL REAGENT	C104-02	6 x 50 mL
CHOLESTEROL REAGENT	C103-01	1 x 250 mL

REAGENT PREPARATION:

The Catachem Cholesterol reagent is packaged ready for use.
No preparation is required.

REAGENT STORAGE AND STABILITY

Store unopened Catachem Cholesterol reagent at 2-8°C (refrigerated).
When stored as directed, the reagent is stable until the expiration date stated on the label.

DIRECTIONS FOR USE

The Catachem Cholesterol Reagent, after opening, is stable on board the analyzer at 2-8°C (refrigerated) for 30 days and capped tightly when not in use.
NOT FOR USE IN UNPROFESSIONAL SETTINGS.

**FOR TECHNICAL ASSISTANCE PLEASE CALL:
(203) 262-0330**

**OR EMAIL:
catachem@catacheminc.com**

www.catacheminc.com