



**URIC ACID REAGENT
C134-0A V134-13**

CONTENTS

PRODUCT NO.

PACKAGING

URIC ACID REAGENT	C134-0A	
URIC ACID ENZYME REAGENT (R1)	C134-04	6 x 24 mL
URIC ACID COLOR REAGENT (R2)	C134-03	6 x 6.5 mL
URIC ACID REAGENT	V134-13	
URIC ACID ENZYME REAGENT (R1)	V134-14	6 x 10 mL
URIC ACID COLOR REAGENT (R2)	V134-15	6 x 3.2 mL

REAGENT PREPARATION

These reagents are packaged ready for use. No preparation is required.

REAGENT STORAGE AND STABILITY

Store unopened reagent at 2-8°C (refrigerated). When stored as directed, the reagent is stable until the expiration date stated on the label.

DIRECTIONS FOR USE

The Catachem Uric Acid Reagents are stable on board the analyzer for 60 days at 2-8°C (refrigerated). Cap reagents when not in use.

FOR TECHNICAL ASSISTANCE PLEASE CALL:

+1 (203) 262-0330 or email catachem@catacheminc.com

Or use the contact form at

www.catacheminc.com