



ASPARTATE AMINOTRANSFERASE (SGOT) REAGENT
V154-0A, C154-0A

<u>CONTENTS</u>	<u>PRODUCT NO.</u>	<u>PACKAGING</u>
AST KIT	V154-0A	
AST SUBSTRATE (R1)	V154-09	3 x 13 mL
AST ACTIVATOR (R2)	V154-10	3 x 2.5 mL

REAGENT PREPARATION

The AST SUBSTRATE (R1) and the AST ACTIVATOR (R2) are packaged ready for use. No preparation is required.

REAGENT STORAGE AND STABILITY

Store unopened reagents at 2-8°C (refrigerated). When stored as directed, the reagents are stable until the expiry date stated on the label.

DIRECTIONS FOR USE

The Catachem AST Test method requires two reagents: AST Substrate as R1 and AST Activator as R2. After opening, these reagents are stable on board the analyzer for 30 days at 2-8°C (refrigerated) for 30 days. Cap reagents tightly when not in use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

**FOR TECHNICAL ASSISTANCE PLEASE CALL:
(203) 262-0330**

**OR EMAIL:
catachem@catacheminc.com**

www.catacheminc.com