



CARBON DIOXIDE (CO₂) REAGENT V206-12

CONTENTS

CARBON DIOXIDE (CO₂) REAGENT
(Cycling)

PRODUCT NO.

V206-12

PACKAGING

6 x 13 mL

REAGENT PREPARATION

Catachem Carbon Dioxide (cycling) is a liquid, ready to use reagent. No preparation is required.

REAGENT STORAGE AND STABILITY

Store unopened reagent at 2-8°C (refrigerated). When stored as directed, the reagent is stable until the expiration date stated on the label.

DIRECTIONS FOR USE

The Catachem Carbon Dioxide Reagent, after opening, is stable on board the analyzer for 15 days at 2-8°C (refrigerated) and tightly capped when not in use.
NOT FOR USE IN UNPROFESSIONAL SETTINGS.

FOR TECHNICAL ASSISTANCE PLEASE CALL:

(203) 262-0330

OR EMAIL:

catachem@catacheminc.com

WWW.CATACHEMINC.COM