



BLOOD UREA NITROGEN (BUN) REAGENT
V264-12, C264-03, C264-04

CONTENTS

PRODUCT NO.

PACKAGING

BLOOD UREA NITROGEN (BUN)	V264-12	6 x 13 mL
BLOOD UREA NITROGEN (BUN)	C264-03	6 x 25 mL
BLOOD UREA NITROGEN (BUN)	C264-04	6 x 50 mL

REAGENT PREPARATION

This reagent is packaged ready for use. No preparation is required.

REAGENT STABILITY AND STORAGE

Store unopened reagent at 2-8°C (refrigerated). When stored as directed, the reagent is stable until the expiration date stated on the label.

DIRECTIONS FOR USE

The Catachem BUN reagent, after opening, is stable on board the analyzer for 30 days at 2-8°C (refrigerated). Cap reagent when not in use.
NOT FOR USE IN UNPROFESSIONAL SETTINGS.

FOR TECHNICAL ASSISTANCE PLEASE CALL:

(203) 262-0330

OR EMAIL:

catachem@catacheminc.com

www.catacheminc.com