



<u>CONTENTS</u>	<u>PRODUCT NO.</u>	<u>PACKAGING</u>
INORGANIC PHOSPHORUS REAGENT	V274-12	6 x 13 mL
INORGANIC PHOSPHORUS BULK	C274-01	500mL

### REAGENT PREPARATION

This reagent is packaged ready for use. No preparation is required.

### REAGENT STORAGE AND STABILITY

Store unopened reagent at 15 - 30°C (room temperature) or at 2 - 8°C (refrigerated).  
When stored as directed, the reagent is stable until the expiration date stated on the label.

### DIRECTIONS FOR USE

The Catachem Inorganic Phosphorus Reagent, after opening, is stable until the expiration date stated on the label. Opened reagent may also be stored on board an automated analyzer at 2 - 8°C (refrigerated) or at 8 - 15°C. Cap reagent when not in use.  
NOT FOR USE IN UNPROFESSIONAL SETTINGS.

## **FOR TECHNICAL ASSISTANCE PLEASE CALL:**

**(203) 262-0330**

**OR EMAIL:**

**[catachem@catacheminc.com](mailto:catachem@catacheminc.com)**

**[www.catacheminc.com](http://www.catacheminc.com)**