



<u>CONTENTS</u>	<u>PRODUCT NO.</u>	<u>PACKAGING</u>
AMYLASE REAGENT KIT	V334-12	6 x 13 mL
AMYLASE REAGENT	C332-01	1 x 250 mL
AMYLASE REAGENT KIT	C334-01	6 x 25 mL

REAGENT PREPARATION

This reagent is packaged ready for use. No preparation is required.

REAGENT STORAGE AND STABILITY

Store unopened reagent at 2 - 8°C (refrigerated). When stored as directed, the reagent is stable until the expiration date stated on the label.

DIRECTIONS FOR USE

The Catachem Amylase Reagent, after opening, is stable on board the analyzer at 2 - 8°C (refrigerated) for 30 days. Cap reagent when not in use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

**FOR TECHNICAL ASSISTANCE PLEASE CALL:
(203) 262-0330**

**OR EMAIL:
catachem@catacheminc.com**

www.catacheminc.com