



**CONTENTS**

**PRODUCT NO.**

**PACKAGING**

LACTATE REAGENT

V452-12

3 x 13 mL

LACTATE REAGENT

C454-01

1 x 25 mL

LACTATE REAGENT KIT A

C452-0A

LACTATE REAGENT

C454-01

6 x 25 mL

**REAGENT PREPARATION**

This reagent is packaged ready for use. No preparation is required.

**REAGENT STORAGE AND STABILITY**

Store unopened reagent at 2-8°C (refrigerated). When stored as directed, reagent is stable until the expiration date stated on the label.

**DIRECTIONS FOR USE**

The Catachem Lactate Reagent, after opening, is stable on board the analyzer at 2-8°C (refrigerated) for 60 days. Keep tightly capped when not in use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

**FOR TECHNICAL ASSISTANCE PLEASE CALL:**

**(203) 262-0330**

**Or use the contact form at**

**[www.catacheminc.com](http://www.catacheminc.com)**