



CALCIUM REAGENT
C295-50, C295-20, C295-05

CONTENTS

PRODUCT NO.

PACKAGING

CALCIUM REAGENT

C295-50

4 x 250 mL

C295-51

1 x 250 mL

CALCIUM REAGENT

C295-20

4 x 50 mL

C295-21

1 x 50 mL

CALCIUM REAGENT

C295-05

2 x 250 mL

C295-51

1 x 250 mL

REAGENT PREPARATION

This reagent is packaged ready for use. No preparation is required.

REAGENT STORAGE AND STABILITY

Store unopened reagent at 2-30 °C (room temperature or refrigerated). When stored as directed, the reagent is stable until the expiration date stated on the label.

DIRECTIONS FOR USE

The Catachem Calcium Reagent, after opening, is stable for 60 days. Opened reagent may also be stored on board an automated analyzer at 2-8 °C (refrigerated) or off an automated analyzer at 8-30 °C for 60 days. Cap reagent when not in use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

FOR TECHNICAL ASSISTANCE PLEASE CALL:

(203) 262-0330

Or use the contact form at

www.catacheminc.com