

**CONTENTS****PRODUCT NO.****PACKAGING****CALCIUM REAGENT****C295-50****4 x 250 mL****C295-51****1 x 250 mL****CALCIUM REAGENT****C295-20****4 x 50 mL****C295-21****1 x 50 mL****CALCIUM REAGENT****C295-05****2 x 250 mL****C295-51****1 x 250 mL****REAGENT PREPARATION**

This reagent is packaged ready for use. No preparation is required.

**REAGENT STORAGE AND STABILITY**

Store unopened reagent at 2-30 °C (room temperature or refrigerated). When stored as directed, the reagent is stable until the expiration date stated on the label.

**DIRECTIONS FOR USE**

The Catachem Calcium Reagent, after opening, is stable for 60 days. Opened reagent may also be stored on board an automated analyzer at 2-8 °C (refrigerated) or off an automated analyzer at 8-30 °C for 60 days. Cap reagent when not in use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

**FOR TECHNICAL ASSISTANCE PLEASE CALL:**

**(203) 262-0330**

**Or use the contact form at**

**[www.catacheminc.com](http://www.catacheminc.com)**