



γ-GLUTAMYL TRANSPEPTIDASE (GGT) REAGENT KIT
V714-0A, C714-0A

<u>CONTENTS</u>	<u>PRODUCT NO.</u>	<u>PACKAGING</u>
GGT (CARBOXY GNA) KIT	V714-0A	
GGT DILUENT (R1)	V714-12	3 x 13 mL
GGT SUBSTRATE (R2)	V714-13	3 x 3.5 mL
GGT (CARBOXY GNA) KIT	C714-0A	
GGT DILUENT (R1)	C714-01	6 x 22.5 mL
GGT SUBSTRATE (R2)	C714-02	6 x 5.5 mL

REAGENT PREPARATION

The Catachem GGT reagents are liquid, ready to use.
No preparation is required.

REAGENT STORAGE AND STABILITY

Store unopened reagents at 2-8°C (refrigerated). When stored as directed, the reagents are stable until the expiration date stated on the label.

DIRECTIONS FOR USE

The Catachem GGT Test method requires two reagents: GGT Diluent as R1 and GGT Substrate as R2. After opening, these reagents are stable on board the analyzer at 2-8°C (refrigerated) for 30 days and tightly capped when not in use.
NOT FOR USE IN UNPROFESSIONAL SETTINGS.

FOR TECHNICAL ASSISTANCE PLEASE CONTACT

Catachem at +1203-262-0330 or Catachem@catacheminc.com