



**CHLORIDE REAGENT KIT
V194-0A, C195-0A, C195-0B**

<u>CONTENTS</u>	<u>PRODUCT NO.</u>	<u>PACKAGING</u>
CHLORIDE REAGENT KIT	V194-0A	
CHLORIDE SAMPLE DILUENT (R1)	V194-01	3 x 13 mL
CHLORIDE COLOR (R2)	V194-02	3 x 4.5 mL
CHLORIDE REAGENT KIT	C195-0A	
CHLORIDE SAMPLE DILUENT (R1)	C195-01	3 x 60 mL
CHLORIDE COLOR (R2)	C195-02	3 x 20 mL
CHLORIDE REAGENT KIT	C195-0B	
CHLORIDE SAMPLE DILUENT (R1)	C195-08	2 x 25 mL
CHLORIDE COLOR (R2)	C195-09	2 x 9.0 mL

REAGENT PREPARATION

The Catachem reagents are packaged ready for use. No preparation of either reagent is required.

REAGENT STORAGE AND STABILITY

Store the unopened reagent kit at 15 - 30°C (room temperature). When stored as directed, both unopened reagents are stable until the expiration date stated on the label.

DIRECTIONS FOR USE

The Catachem Chloride test method requires two reagents: the Buffer (V194-01) as R1 and the Color (V194-02) as R2. After opening, these reagents are stable at 2 - 30°C for at least 30 days. Cap when not in use. NOT FOR USE IN UNPROFESSIONAL SETTINGS.

**FOR TECHNICAL ASSISTANCE PLEASE CALL:
(203) 262-0330**

**OR EMAIL:
catachem@catacheminc.com
www.catacheminc.com**